
RI EZ-Squeeze™

Instructions For Use

CE R_X Only
2797

+45 46 79 02 02 | sales@coopersurgical.com | fertility.coopersurgical.com

Research Instruments Ltd,
Bickland Industrial Park, Falmouth, Cornwall TR11 4TA, UK
CooperSurgical Distribution B.V. Celsiusweg 35, 5928 PR
Venlo, The Netherlands

DEVICE DESCRIPTION

EZ-Squeeze is a plastic pipette tip fitted with a silicone bulb. The EZ-Squeeze is available in a range of internal diameter sizes. Sizes 200µm, 290µm and 600µm are suitable for specimen handling. Sizes 131µm, 135µm, 145µm, 155µm and 170µm are suitable for oocyte denudation.

Denudation may be performed on oocytes in preparation for Intracytoplasmic Sperm Injection (ICSI) and In Vitro Fertilisation (IVF) techniques. It allows easier access of the sperm to the zona pellucida and better visibility, and therefore more effective manipulation of the injection pipette whilst conducting ICSI procedures. The technique of denuding is where the oocyte is repeatedly aspirated and expelled through pipette tips of incrementally smaller internal diameters to remove the outer layers of cumulus and corona cells. Generally two sizes of pipette are used to achieve this. The sizes selected depend on the size of the oocyte, and user preference and practice.

CAUTION

- **Caution:** US Federal Law restricts this device to sale by or on the order of a physician.
- Do not re-use - reuse poses a high risk of cross contamination.

Caution: These tools are intended for one-time use only and must remain covered until ready to use in order to maintain both sterility and prevent contamination. The user is warned to follow safe, established laboratory practices in order to prevent accidental DNA contamination after the EZ-Squeeze has been removed from its protective covering.

2

- EZ-Squeeze should only be used by qualified personnel who have been approved by the relevant authorised bodies.
- The introduction of EZ-Squeeze for routine clinical use should be risk assessed prior to use.

Caution: This tool may be used to handle polar body, blastomere or trophectoderm biopsy samples for pre-implantation genetic diagnosis (PGD) on the genetic material in the biopsied cell(s). Test for PGD are currently developed and their performance characteristics are determined by individual laboratories for their own use. The performance of these tests may vary depending on the particular assay and disease evaluated. Currently, these tests have not been cleared or approved by the US Food and Drug Administration (FDA).

This product has a 5 year shelf life.

INTENDED USE

For denudation i.e. to remove the cumulus from an oocyte prior to Intracytoplasmic Sperm Injection (ICSI) and In Vitro Fertilization (IVF) and for handling gametes, embryos and biopsied cells (polar bodies, blastomeres and trophectoderm) during Assisted Reproductive Technology (ART) procedures.

CONTRAINDICATIONS

EZ-Range pipettes are not intended for biopsy of cells from oocytes or embryos.

3

INTENDED PATIENT POPULATION

For those who require assisted reproductive treatments to improve their chances of a healthy pregnancy and outcome.

MEDICAL CONDITION

The medical condition is infertility and genetic disease.

CLINICAL BENEFIT

The clinical benefit of using this medical device during ART procedures is to improve the patient's chances of a viable pregnancy and long-term health of the offspring.

INTENDED ENVIRONMENT OF USE

For use in a controlled clinical environment where Assisted Reproductive Techniques (ART) are used.

INTENDED USER

Trained clinical professional working in Assisted Reproductive Techniques (ART) with locally relevant qualifications.

DISPOSAL

- Remove EZ-Squeeze from the handle.
- Dispose of the EZ-Squeeze and packaging in a suitable container.

4

BEFORE USE

- Check that the packaging has not been damaged. If the sterile barrier has been compromised, discard.
- Check expiry date on pouch. If out of date, discard.
- Select the appropriate pipette tip size to ensure minimal distortion and damage to the sample i.e. the internal diameter of the pipette tip should be appropriate for the size of the sample to be denuded or transferred. This can be done by placing the pipette tip next to the sample to determine suitability.
- Check the pipette tip under your microscope to ensure it is clear. Rotate the tip 90° and check again. If it is not clear, discard.
- To avoid creating bubbles in the media first apply pressure to the silicone bulb before immersing the tip.
- Rinse the tip by depressing the membrane and immersing the tip in some media. Then slowly release the pressure to aspirate media into the pipette tip. Depress the membrane again to expel the media. Repeat this action a few times to complete the rinsing process.
- If using for the first time (before attempting clinical use), familiarise yourself with the device using plain water or media, and discarded eggs or embryos. Dispose of the practice tip after use.

3. Bend back the tray to expose the bulb end of the EZ-Squeeze and remove.
4. Insert the silicone bulb (fitted with the pipette tip) into a compatible handle e.g. EZ-Squeeze handle (See Fig.2)


Fig.1


Fig.2

HOW TO USE

1. For denudation, transfer the sample into media with an enzyme constituent, hyaluronidase, which promotes the denuding process.

5

2. Apply pressure to the bulb at the exposed section of the handle with finger or thumb pressure and immerse the tip into the media. If working under oil, it is important to apply positive pressure when passing through the oil layer.
3. Once the sample is located, slowly release the pressure on the bulb to aspirate it with some media into the EZ-Squeeze. Transfer the sample at this stage if required.
4. Apply pressure to the bulb again to expel the sample /media.
5. Repeat this action a sufficient number of times to denude the sample.
6. Always keep the media level visible in the EZ-Squeeze.

OBLIGATION TO INFORM

Any serious incident that has occurred in relation to this device should be reported to CooperSurgical via phone number +1 203-601-5200 Ext 3100 or by email at ProductSurveillance@coopersurgical.com and to the local Health Authority in your country. A serious incident may have caused or contributed to a death, a delay in a procedure which resulted in death or serious injury, or a malfunction that could have caused an adverse event.

6

QUALITY ASSURANCE

- Sterilised by gamma irradiation to a Sterility Assurance Level (SAL) of 10⁻⁶.
- Non-embryotoxic, Pass Level: 1-Cell MEA: ≥80% hatched blastocysts at 120h.
- Non-Pyrogenic (Limulus Amebocyte Lysate [LAL] Test), Pass Level: <20 EU/device.
- Certificate of analysis available on request.
- Depending on the size of the tip the volumetric capacity of the EZ-Squeeze is 55.6µl - 60.6µl.

SYMBOL GLOSSARY Source: ISO 15223-1


	Consult instructions for use		Do not use if packaging is damaged
	Do not re-use		Sterilised using gamma irradiation
	Keep away from direct sunlight (UV)		80% 20% Store between 20% and 80% relative humidity limits
	Batch code		30°C 10°C Store between 10°C and 30°C temperature limits
	Use-by-date		Date of manufacture
	Manufacturer		Unique Device Identifier
	Medical Device		Authorised representative in the European Community
	Catalogue or Part number		

PREPARATION

1. Remove one blister pack tray from the box and snap off a single EZ-Squeeze unit.
2. Peel back the lid half way, starting at the wide end (See Fig.1).

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