**EC REP** CopperSurgical Distribution B.V.  
Celsiusweg 35, 5928 PR Venlo, The Netherlands  
LifeGlobal Group, LLC, 393 Soundview Rd, Guilford, CT 06437 US  
T: 1800-648-1151 Int'l: +45 469200  
sales@coopersurgical.com Orders@origo.us.com

## Instructions for the Use of HTF w/ HEPES

(Catalogue Numbers: GMHH-100, GMHH-250, GMHH-500)

### PRECAUTIONS AND WARNINGS

- Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
- Caution: The user should read and understand the Instructions for Use, Precautions and Warnings, and be trained in the correct procedure before using the HTF w/ HEPES for retrieval of human oocytes.
- Not to be used for injection.
- Do not re-sterilize.
- Do not use the product if:
  - the product packaging appears damaged or if the seal is broken
  - the product becomes discolored, cloudy, or shows evidence of particulate matter
- HTF w/ HEPES contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- HTF w/ HEPES contains only a low concentration of sodium bicarbonate and should not be gassed, or used for culture under CO<sub>2</sub>.
- To avoid problems with contamination, practice aseptic techniques.
- Discard unused medium within 7 days of opening. Do not use after expiry date.

### GENERAL INFORMATION

#### Indications for Use

For handling, washing & manipulation of embryos & oocytes.

#### Storage and Shelf Life

Store at 2-8°C and protected from light. Fourteen (14) weeks from the date of manufacture.

**Composition**

A HEPES-buffered medium is required when gametes and embryos are handled outside of a CO<sub>2</sub> incubator.

Sodium Chloride, Potassium Chloride, Calcium Chloride, Potassium Phosphate, Magnesium Sulfate, Sodium Bicarbonate, Glucose, Lactate Na Salt, Sodium Pyruvate, Phenol Red, HEPES, Gentamicin Sulfate\* (10 µg/ml)

\*from therapeutic-grade source material

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### INSTRUCTIONS FOR USE

#### Media Required for Oocyte Retrieval

- Needle-Wash Medium: Use HTF w/ HEPES to wash the collection needles prior to oocyte retrieval.
- Follicle-Flush Medium: Use HTF w/ HEPES to flush the follicles, if appropriate.
- Oocyte Wash Medium: Use HTF w/ HEPES supplemented with protein to wash the oocytes after retrieval.
- Oocyte Culture Medium: Use global® for fertilization supplemented with protein for the final wash and culture of the oocytes prior to fertilization.

The procedures described below have been found to be effective for oocyte retrieval and washing, holding during ICSI procedures, and for embryo holding during assisted hatching at Day 3. Every laboratory must define and optimize its own procedures.

The day before oocyte retrieval, aliquot appropriate amount of oil, into a 14 ml Falcon tube, and place it in a warming block (36-38°C) to equilibrate overnight.

Prepare oocyte-wash medium by supplementing the medium with either Human Serum Albumin (HSA) or LifeGlobal® Protein Supplement to achieve desired % (v/v) of protein supplementation.

- Approximately 2 hours prior to oocyte retrieval, use a thermometer to confirm that the tube warmer is at 36-38°C.
- Place a tube of oocyte-wash for each patient into the warming block.
- Place 10-12.100 mm Falcon dishes on the warming surface of the laminar flow hood.
- Prepare dishes containing oocyte culture medium for final oocyte wash and culture.
- Just prior to oocyte retrieval, prepare the initial oocyte wash dish by pouring of oocyte-wash medium into a warmed 60 mm Falcon dish. Cover the oocyte-wash medium with approximately 5.5 ml of pre-equilibrated oil and place the dish on the warming surface of the laminar flow hood.
- During oocyte retrieval, search for the oocytes under a dissecting microscope. Using a sterile pipette, transfer each oocyte into the initial wash dish containing oocyte-wash medium.
- Dissect each oocyte using two 25 g 11/2" needles to remove the degenerated cumulus cells, debris and blood.
- Transfer all oocytes to a dish of the oocyte culture medium for final washing. Pipette all the oocytes a few times to remove traces of oocyte-wash medium.
- Transfer each dissected oocyte into fresh culture media, according to your standard laboratory practice, and place in the incubator until fertilization or ICSI.

### SYMBOLS

STERILE A	RX Only	REF	LOT	■ i	■	■ ☼
Stérile Usine Aspéticas Techniques	By Prescription Only	Catalogue Number	Batch Code	Consult Instructions For Use	Manufacturer	Keep Away from Sunlight
0-8°C	EC REP					
Temperature Limitation	Authorized Representative in the European Community	Use By	GSI DataMatrix Barcode	Do Not Resterilize	Do Not Use If Package is Damaged	

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