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
Quinn's Advantage™ Sperm Freezing-Medium

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Gewicht: 23 kg (netto) 26 kg (brutto)

Literatur

1. DeRosa N, Pooley K, Kohut T, Dissing M, Campbell B, Kirkman-Brown J. Synergistic role of bicarbonate and pH on sperm motility and velocity in sperm preparations. *Birmingham: British Fertility Society Association of Clinical Embryologists Society for Reproduction and Fertility*. 2015 S. 70.
2. Tomlinson M. Optimizing Therapeutic Sperm Washing-Medium: Why are there clear differences in sperm progression and velocity between products? *CooperSurgical ART Scientific*, Ausgabe v. 6. April 2020.
3. Achikanu C, Pendekanti V, Teague R, Publicover S. Effects of pH manipulation, CatSper stimulation and Ca²⁺-store mobilization on [Ca²⁺]_i and behaviour of human sperm. *Hum Reprod*. 01. Okt. 2018; 33(10): 1802-1811.
4. Abou-haila A, Tulsiani DR. Signal transduction pathways that regulate sperm capacitation and the acrosome reaction. *Arch Biochem Biophys*. 01. Mai 2009; 485(1):72-81.
5. Orsi NM, Leese HJ. Protection against reactive oxygen species during mouse preimplantation embryo development: role of EDTA, oxygen tension, catalase, superoxide dismutase and pyruvate. *Mol Reprod Dev*. Mai 2001; 59(1): 44-53.
6. Wu X, Dai H, Liu L, Xu C, Yin Y, Yi J, Bielec MD, Han Y, Li S. Citrate reduced oxidative damage in stem cells by regulating cellular redox signaling pathways and represent a potential treatment for oxidative stress-induced diseases. *Redox Biol*. Febr. 2019; 21; 101057
7. Yun JI, Gong SP, Song YH, Lee ST. Effects of combined antioxidant supplementation on human sperm motility and morphology during sperm manipulation in vitro. *Fertil Steril*. Aug. 2013;100(2):373-8
8. Holmes E, Björndahl L, Kvist U. Hypotonic challenge reduces human sperm motility through coiling and folding of the tail. *Andrologia*. Dez. 2020; 52(11): e13859
9. Parmegiani L, Cognigni GE, Ciampaglia W, Pocognoli P, Marchi F, Filicori M. Efficiency of hyaluronic acid (HA) sperm selection. *J Assist Reprod Genet*. Jan. 2010; 27(1): 13-6.
10. Parmegiani L, Cognigni GE, Bernardi S, Troilo E, Ciampaglia W, Filicori M. "Physiologic ICSI": hyaluronic acid (HA) favors selection of spermatozoa without DNA fragmentation and with normal nucleus, resulting in improvement of embryo quality. *Fertil Steril*. Febr. 2010; 93(2): 598-604
11. Huszar G, Jakab A, Sakkas D, Ozenci C., Cayli S., Delpiano E., Ozkavukcu S. Fertility testing and ICSI sperm selection by hyaluronic acid binding: clinical and genetic aspects. *Reprod BioMed Online* 2007; 14(5): 650-663.
12. West R, Coomarasamy A, Frew L, Hutton R, Kirkman-Brown J, Lawlor M, Lewis S, Partanen R, Payne-Dwyer A, Román-Montañana C, Torabi F, Tsagdi S, Miller D. Sperm selection with hyaluronic acid improved live birth outcomes among older couples and was connected to sperm DNA quality, potentially affecting all treatment outcomes. *Hum Reprod*. 30. Mai 2022; 37(6): 1106-1125.
13. WorriLOW KC, Eid S, Woodhouse D, Perloe M, Smith S, Witmyer J, Ivani K, Khoury C, Ball GD, Elliot T, Lieberman J. Use of hyaluronan in the selection of sperm for intracytoplasmic sperm injection (ICSI): significant improvement in clinical outcomes--multicenter, double-blinded, and randomized controlled trial. *Hum Reprod*. Febr. 2013; 28(2): 306-14.
14. Miller D, Pavitt S, Sharma V, Forbes G, Hooper R, Bhattacharya S, Kirkman-Brown J, Coomarasamy A, Lewis S, Cutting R, Brison D, Pacey A, West R, Brian K, Griffin D, Khalaf Y. Physiological, hyaluronan-selected intracytoplasmic sperm injection for infertility treatment (HABSelect): a parallel, two-group, randomised trial. *Lancet*. 2. Febr. 2019; 393 (10170): 416-422.
15. Scaruffi P, Bovis F, Casciano I, Maccarini E, De Leo C, Gazzo I, Massarotti C, Sozzi F, Stigliani S, Anserini P. Hyaluronic acid-sperm selection significantly improves the clinical outcome of couples with previous ICSI cycles failure. *Andrology*. Mai 2022; 10(4): 677-685.
16. Fen C.T.C., Lee S. N, Lim M. N, Yu S. L. Relationship between Sperm Hyaluronan-Binding Assay (HBA) Scores on Embryo Development, Fertilisation, and Pregnancy Rate in Patients Undergoing Intra-Cytoplasmic Sperm Injection (ICSI). *Proceedings of Singapore Healthcare*. 2013.22.N2: 120-124.
17. Nijs M, Creemers E, Cox A, Janssen M, Vanheusden E, Van der Elst J, Ombelet W. Relationship between hyaluronic acid binding assay and outcome in ART: a pilot study. *Andrologia*. Okt. 2010; 42(5): 291-6.
18. Nazif MS, Rehman ZU, Khan H, Khan FA, Hussain T, Ahmad A, Farmanullah, Husnain A, Muhammad S, Murtaza G, Gang L. Glycine Improved Cryopreserved Spermatozoa Quality in Achai Bull. *Biomed Res Int*, 4. Aug. 2022; 2022: 8282387
19. Renard P, Grizard G, Grieveau JF, Sion B, Boucher D, Le Lannou D. Improvement of motility and fertilization potential of postthaw human sperm using glutamine. *Cryobiology*. Juni 1996; 33(3): 311-9.
20. Hezavehei M, Sharafi M, Kouchesfahani HM, Henkel R, Agarwal A, Esmaeili V, Shahverdi A. Sperm cryopreservation: A review on current molecular cryobiology and advanced approaches. *Reprod Biomed Online*. Sept. 2018; 37(3): 327-339.
21. Kuo YL, Tzeng WL, Chiang HK, Ni RF, Lee TC, Young ST. New system for long-term monitoring of sperm motility: EDTA effect on semen. *Arch Androl*. Sept./Okt. 1998-; 41(2): 127-33

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