English

SAGE Media™ PureCeption™ Sperm Separation Media

For laboratory procedures only other uses must be qualified by the end user.

REF	Product Description	Number	Unit Size
PureCeption™	ART-2024	12 x 12 mL Bi-Layer Kit	
24-Determination			

ART-2024 is assembled with 6 x 12 mL PureCeption™ 40% (v/v) and 6 x 12 mL PureCeption™ 80% (v/v).

PureCeption™ is assembled in a 24-determination kit. Each determination uses 2.0 mL of Upper Phase (40%) and 2.0 mL of Lower Phase (80%).

MATERIALS PROVIDED IN THE SPERM SEPARATION MEDIA KIT

- Lower Phase Gradient: 80% PureCeption™ with HEPES-buffered Human Tubal Fluid (HTF-HEPES).
- Upper Phase Gradient: 40% PureCeption™ with HEPES-buffered Human Tubal Fluid (HTF-HEPES).

PURECEPTION™ GRADIENTS ARE ALSO AVAILABLE IN 100 mL VOLUMES

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Gradients containing different percentages of PureCeption™ Sperm Separation Media can be made by diluting the PureCeption™ 100% Isotonic Solution with appropriate proportions of Quinn's Sperm Washing Medium (ART-1005/1006).

Each laboratory should make its own determination for which concentration of the PureCeption™ Sperm Separation Medium they utilize for the Upper or Lower Phase based on the quality of the semen sample.

INTENDED USE

These products are intended to facilitate highly efficient separation of motile sperm from the ejaculate of most semen specimens.

PRODUCT DESCRIPTION

A 40% (v/v), 80% (v/v), or 100% (v/v) sterile colloid suspension of silica particles stabilized with covalently bound hydrophilic silane supplied with HEPES-buffered Human Tubal Fluid (HTF). The components of this kit will allow for the highly efficient separation of motile sperm from the ejaculate of most semen specimens. It is fast, cost-effective, and needs very little equipment or expertise to perform.

This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

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PRECAUTIONS AND WARNINGS

PureCeption™ comes packaged with tamper-proof seals and caps. If the seal is broken or the cap loose, do not use the product.

PureCeption™ has a naturally opalescent appearance. Do not use product if it shows any evidence of particulate matter or contamination. This may be evident by extreme cloudiness or discoloration.

PureCeption™ should remain tightly capped when placed in a CO₂ incubator to avoid pH changes.

Single Use: To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and embryo quality during *in-vitro* culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

Note: Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

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1. ART-1005/1006 Quinn's Sperm Washing Medium
2. Sterile disposable polystyrene 15 mL conical centrifuge tube.
3. Sterile 5 mL disposable pipettes (graduated/transfer) or 3-cc syringes with 1.5/21g needles.
4. Centrifuge/fixed or horizontal. Must be able to operate for up to 30 minutes at 250 x g to 750 x g.
5. 37 °C incubator or water bath.
6. Counting chamber.
7. Microscope with x10 and x20 objectives.

QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

QUALITY CONTROL

All PureCeption™ solutions are tested using one-cell mouse embryo culture and endotoxin assays to ensure quality and safety. However, it is recommended that inhouse quality control be performed with each lot.

USER QUALITY CONTROL

1. Examine the PureCeption™ solutions. Do not use any medium that appears discolored. The PureCeption™ Upper Phase (40%), Lower Phase (80%), and Isotonic Solution (100%) should have a normal opalescent appearance.
2. To avoid contamination:
 a. Enter all bottles using sterile technique within a suitable sterile environment.
 b. Do not use the same sterile pipette or syringe to reenter a bottle of medium.
 c. When entering a bottle through the stopper via a needle, use a sterile needle.
 d. Wipe the stopper with alcohol and air dry. Use a new sterile needle for each gradient.

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NOTE: Use the PureCeption™ gradient within 1 hour after creating it. Otherwise, the two phases will eventually blend into each other and a sharp interface will not exist. Considerable particulate debris accumulates at the Upper/Lower Phase interface. This is considered a valuable part of the PureCeption™ system.

3. Gently place up to 2.5 mL of liquefied semen onto the Upper Phase using a pipette or syringe.
4. Centrifuge for 20 minutes at 350 x g to 400 x g (or up to 750 x g for highly viscous samples). When first centrifugation is completed, you may not be able to see a pellet but it is essential to continue the procedure.
5. Following centrifugation, remove all layers except the lowest portion (approximately 0.3 mL).
6. Transfer the pellet and remaining suspension liquid into a sterile conical tube with a sterile disposable pipette. This will avoid resuspending the pellet with any debris from the wall of the conical tube from the initial centrifugation.
7. Add 2 to 3 mL of Quinn's Sperm Washing Medium and resuspend the pellet.

DIRECTIONS FOR USE FOR PROCEDURE FOR FRESH SEMEN SAMPLES

1. Bring all appropriate PureCeption™ sperm media and the semen sample to room temperature before using. This will avoid a cold "shock" to the spermatozoa.

2. To prepare the single layered gradient:
 a. Transfer 2.0 mL of Lower Phase (PureCeption™ 80%) into the conical tube.

b. Using a pipette or syringe, layer 2.0 mL of Upper Phase (PureCeption™ 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.

c. Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly disperse the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single layered gradient is stable for 1 hour.

8. Centrifuge for 4 to 8 minutes at 250 x g. The higher number of sperm will require the maximum 8-minute centrifugation time to ensure a complete and thorough sperm wash.

9. Remove supernatant and resuspend the pellet in 0.5 mL of Suitable Washing Solution.

10. Sample is now ready for recount and estimate of motility.

11. Sample is now ready for recount and estimate of motility.

MATERIALS REQUIRED BUT NOT INCLUDED

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