## English

REF	Number	Unit Size
PureCepcion <sup>™</sup> 40%	ART-2040	100 mL
PureCepcion <sup>™</sup> 80%	ART-2080	100 mL
PureCepcion <sup>™</sup> 100%	ART-2100	100 mL

For laboratory procedures only; other uses must be qualified by the end user.

REF	Number	Unit Size
PureCepcion <sup>™</sup> 40%	ART-2024	12 x 12 mL

ART-2024 is assembled with 6 x 12 mL PureCepcion<sup>™</sup> 40% Upper Phase and 6 x 12 mL PureCepcion<sup>™</sup> 80% Lower Phase.

PureCepcion<sup>™</sup> is assembled in a 24-Determination kit. Each determination uses 2.0 mL of Upper Phase (40%) and 2.0 mL of Lower Phase (80%).

**MATERIALS PROVIDED IN THE SPERM SEPARATION MEDIA KIT**

- Lower Phase Gradient: 80% PureCepcion<sup>™</sup> with HEPES-buffered Human Tubal Fluid (HTF-HEPES).
- Upper Phase Gradient: 40% PureCepcion<sup>™</sup> with HEPES-buffered Human Tubal Fluid (HTF-HEPES).

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**PURECEPTION<sup>™</sup> GRADIENTS ARE ALSO AVAILABLE IN 100 mL VOLUMES**

REF	Number	Unit Size
PureCepcion <sup>™</sup> 40%	ART-2040	100 mL
PureCepcion <sup>™</sup> 80%	ART-2080	100 mL
PureCepcion <sup>™</sup> 100%	ART-2100	100 mL

Gradients containing different percentages of PureCepcion<sup>™</sup> can be made by diluting the PureCepcion<sup>™</sup> 100% Isotonic Solution with appropriate proportions of Quinn's Sperm Washing Medium (ART-1005/1006).

**Each laboratory should make its own determination for which concentration of the PureCepcion<sup>™</sup> they utilize for the Upper or Lower Phase based on the quality of the semen sample.**

**INTENDED USE**

These products are intended to facilitate highly efficient separation of motile sperm from the ejaculate of most semen specimens.

**PRODUCT DESCRIPTION**

A 40% (v/v), 80% (v/v), or 100% (v/v) sterile colloid suspension of silica particles stabilized with covalently bound hydrophilic silane supplied with HEPES-buffered Human Tubal Fluid (HTF). The components of this kit will allow for the highly efficient separation of motile sperm from the ejaculate of most semen specimens. It is fast, cost-effective, and needs very little equipment or expertise to perform.

This product contains 10 mg/L of gentamicin, an aminoglycoside antibiotic.

## PRECAUTIONS AND WARNINGS

PureCepcion<sup>™</sup> comes packaged with tamper-proof seals and caps. If the seal is broken or the cap loose, do not use the product. PureCepcion<sup>™</sup> has a naturally opalescent appearance. Do not use product if it shows any evidence of particulate matter or contamination. This may be evidence by extreme cloudiness or discoloration.

PureCepcion<sup>™</sup> should remain tightly capped when placed in a CO<sub>2</sub> incubator to avoid pH changes.

**Single Use:** To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

The products are aseptically processed and supplied sterile.

**Note:** Embryo is considered a general term. More precisely, SAGE<sup>™</sup> considers the period of time initiating when a single diploid cell results from the fusion of male and female genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day

5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

**QUALITY ASSURANCE**

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

**QUALITY CONTROL**

All PureCepcion<sup>™</sup> solutions are tested using one-cell mouse embryo culture and endotoxin assays to ensure quality and safety. However, it is recommended that inhouse quality control be performed with each lot.

**USER QUALITY CONTROL**

- Examine the PureCepcion<sup>™</sup> solutions. Do not use any medium that appears discolored. The PureCepcion<sup>™</sup> Upper Phase (40%), Lower Phase (80%), and Isotonic Solution (100%) should have a normal opalescent appearance.
- To avoid contamination:
  - Enter all bottles using sterile technique within a suitable sterile environment.
  - Do not use the same sterile pipette or needle to enter a bottle of medium.
  - When entering a bottle through the stopper via a needle, use a sterile needle.
  - Wipe the stopper with alcohol and air dry. Use a new sterile needle for each gradient.

## MATERIALS REQUIRED BUT NOT INCLUDED

- ART-1005/1006 Quinn's Sperm Washing Medium
- Sterile disposable polystyrene 15 mL conical centrifuge tube.
- Sterile 5 mL disposable pipettes (graduated/transfer) or 3-cc syringes with 1.5/2/1g needles.
- Centrifugefixed or horizontal. Must be able to operate for up to 30 minutes at 250 x g to 750 x g.
- 37 °C incubator or water bath.
- Counting chamber.
- Microscope with x10 and x20 objectives.

**DIRECTIONS FOR USE FOR PROCEDURE FOR FRESH SEMEN SAMPLES**

- Bring all appropriate PureCepcion<sup>™</sup> sperm media and the semen sample to room temperature before using. This will avoid a cold "shock" to the spermatozoa.
- To prepare the single bilayered gradient:
  - Transfer 2.0 mL of Lower Phase (PureCepcion<sup>™</sup> 80%) into the conical tube.
  - Using a pipette or syringe, layer 2.0 mL of Upper Phase (PureCepcion<sup>™</sup> 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.
  - Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly disperse the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single bilayered gradient is stable for 1 hour.

**NOTE: Use the PureCepcion<sup>™</sup> gradient within 1 hour after creating it. Otherwise, the two phases will eventually blend into each other and a sharp interface will not exist. Considerable particulate debris accumulates at the Upper/Lower Phase interface. This is considered a valuable part of the PureCepcion<sup>™</sup> system.**

- Gently place up to 2.5 mL of liquefied semen onto the Upper Phase using a pipette or syringe.
- To prepare the single bilayered gradient:
  - Transfer 1.5 mL of Lower Phase (PureCepcion<sup>™</sup> 80%) into a conical tube.
  - Using a pipette or syringe, layer 1.5 mL of Upper Phase (PureCepcion<sup>™</sup> 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.
  - Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly disperse the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single bilayered gradient is stable for 1 hour.

- Bring all appropriate PureCepcion<sup>™</sup> sperm media and the semen sample to room temperature before using. This will avoid a cold "shock" to the spermatozoa.
- To prepare the single bilayered gradient:
  - Transfer 2.0 mL of Lower Phase (PureCepcion<sup>™</sup> 80%) into the conical tube.
  - Using a pipette or syringe, layer 2.0 mL of Upper Phase (PureCepcion<sup>™</sup> 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.
  - Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly disperse the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single bilayered gradient is stable for 1 hour.

**PROCEDURE FOR FROZEN SPERM SAMPLES**

Frozen-thawed sperm tend to be less vigorous than fresh sperm. For this reason, use 1.5 mL of Upper and Lower Phases instead of 2.0 mL volumes recommended for fresh sperm.

- Bring all appropriate PureCepcion<sup>™</sup> sperm media and the semen sample to room temperature before using. This will avoid a cold "shock" to the spermatozoa.
- To prepare the single bilayered gradient:
  - Transfer 1.5 mL of Lower Phase (PureCepcion<sup>™</sup> 80%) into a conical tube.
  - Using a pipette or syringe, layer 1.5 mL of Upper Phase (PureCepcion<sup>™</sup> 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.
  - Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly disperse the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single bilayered gradient is stable for 1 hour.

**NOTE: Use the PureCepcion<sup>™</sup> gradient within 1 hour after creating it. Otherwise, the two phases will eventually blend into each other and a sharp interface will not exist. Considerable particulate debris accumulates at the Upper/Lower Phase interface. This is considered a valuable part of the PureCepcion<sup>™</sup> system.**

- Gently place the thawed semen sample onto the Upper Phase using a pipette or syringe.
- Centrifuge for 20 minutes at 350 x to 400 x g. When first centrifugation is completed,

you may not be able to see a pellet but it is essential to continue the procedure.

- Following centrifugation, remove all layers except lowest portion (approximately 0.5 mL).
- Transfer the pellet and remaining suspension liquid into a sterile conical tube with a sterile disposable pipette. This will avoid resuspending the pellet with any debris from the wall of the conical tube from the initial centrifugation.
- Add 2 to 3 mL of Quinn's Sperm Washing Medium and resuspend the pellet.
- Centrifuge for 8 minutes at 250 x g.
- Remove supernatant and resuspend the pellet to 0.5 mL of Suitable Washing Solution.
- Sample is now ready for estimate of count and motility.

**Each laboratory should make its own determination of which medium and protocol to use for each particular procedure.**

Information on specific aspects of IVF and embryo culture is available in our Product Catalog

**TROUBLESHOOTING**

- Occasionally, samples that do not liquefy properly and remain too viscous to pass through the gradient will be encountered. Increasing the centrifugal force up to but no more than 750 x will aid in separating the sperm in these cases. Pelleting motile sperm from highly viscous semen does not usually present a problem using the PureCepcion<sup>™</sup> system.
- The most important semen parameter contributing to a pellet of highly motile sperm is sperm progression. The higher the number

of progressively motile sperm in the initial sample, the greater the number of sperm you will have in the final pellet. If sperm prove motility or progression is below WHO guidelines for normal semen parameters, you may be able to compensate for shortcomings in the initial semen analysis by using 1.0 mL of both 40% and 80% PureCepcion<sup>™</sup> instead of 2.0 mL of each. You may also want to leave the last 0.5 to 0.7 mL of the post-centrifugation Lower Phase above the pellet and wash the sperm caught in transit through this portion of the Lower Phase. If this is done, be sure to double the volume of Suitable Washing Solution in order to dilute out this remaining Lower Phase material.

**STORAGE INSTRUCTIONS AND STABILITY**

All PureCepcion<sup>™</sup> components must be stored at 2 °C to 8 °C. The product is stable until the expiration date shown on the label.

A. Remove desired volume of product using aseptic procedures.

B. Once removed, do not return any volume of product to the original container.

C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

**RELATED PRODUCTS**

ART-1005/1006 Quinn's Sperm Washing Medium

SAGE<sup>™</sup> In Vitro Fertilization<sup>™</sup> has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE<sup>™</sup> Support Line.

PureCepcion<sup>™</sup> is a trademark of CooperSurgical, Inc.

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