EC REP CopperSurgical Distribution B.V.
Celsiusweg 35, 5928 PR Venlo, The Netherlands
LifeGlobal Group, LLC, 393 Soundview Rd, Guilford, CT 06437 US
T: 1800-648-1151 Int'l: +45 46790200
sales@coopersurgical.com Orders@origo.us

EN

Instructions for the Use of HTF

(Catalogue Numbers: GMHT-100, GMHT-250, GMHT-500)

PRECAUTIONS AND WARNINGS

- Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
- Caution:** The user should read and understand the Instructions for Use, Precautions and Warnings, and be trained in the correct procedure before use for the Culture of cleavage-stage human embryos (Day 1-3).
- Not to be used for injection.
- Do not resterilize.
- Do not use the product if:
 - the product packaging appears damaged or if the seal is broken
 - the expiry date has been exceeded
 - the product becomes discolored, cloudy, or shows evidence of particulate matter
- HTF contains antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.
- To avoid problems with contamination, practice aseptic techniques.
- Discard unused medium within 7 days of opening. Do not use after expiry date.

GENERAL INFORMATION

- Indications for Use**
Culture of cleavage-stage human embryos (Day 1 to Day 3).

Storage and Shelf Life

Store at 2-8°C and protected from light. Ten (10) weeks from the date of manufacture.

- Composition**
A bicarbonate-buffered medium replete with glucose, lactate, and pyruvate is required to support the growth and development of human embryos in vitro.

Sodium Chloride, Potassium Chloride, Calcium Chloride, Potassium Phosphate, Magnesium Sulfate, Sodium Bicarbonate, Glucose, Lactate Na Salt, Sodium Pyruvate, Phenol Red, Gentamicin Sulfate* (10 µg/ml)

*from therapeutic-grade source material

QUALITY CONTROL SPECIFICATIONS

Assay (performed for each batch)	Specification
Physicochemical Tests	
pH (with 5% CO ₂)	7.2-7.4
Osmolalität	280-292 mOsm
Biological Tests	
Endotoxin (LAL)	≤ 0.5 EU/ml
Test de stérilité (dépistage bactérien et fongique, SAL 10 ⁻³)	RÉUSSI
Biological Assays	
1-cell Mouse Embryo Assay (% expanded blastocysts at 96 h of culture)	≥ 80%

Special Note on the CO₂ Concentration in the Incubator: In most cases, a 5-7% concentration of CO₂ in the incubator will produce a pH of 7.2 to 7.4 in HTF. However, the exact concentration of CO₂ required to produce the optimum pH of approximately 7.30 (7.27-7.33) depends on several factors, including the physical characteristics of incubator and the altitude. Consequently, we strongly recommend that each laboratory determine and use the concentration of CO₂ that is required to produce a pH of 7.30 in HTF.

INSTRUCTIONS FOR USE

The procedures described below have been found to be effective for the culture of human embryos. Every laboratory must define and optimize its own procedures.

After each time the original bottle is opened recap the bottle tightly and store at 2-8°C, protected from light.

Twenty-four (24) hours prior to the use of HTF for embryo culture, supplement the medium with either Human Serum Albumin (HSA) or LifeGlobal® Protein Supplement to achieve desired % (v/v) of protein supplementation.

CO₂ and Temperature Equilibration

Method 1: The day before embryo culture:
1. Slowly gas the supplemented HTF in the 50 ml tissue culture flask for approximately 15 seconds with 5:5:90 Blood Gas Mixture using a sterile Pasteur pipette.

2. Allow the bubbles to rise to the neck of the flask.

3. Cap the culture flask tightly and store in the refrigerator at 2-8°C. Over night the gas bubbles will disperse into the media, bringing the medium to the appropriate pH (7.2-7.4).

4. The next day, 1 hour prior to embryo culture, prepare microdrops of the gassed, supplemented HTF and place in the incubator for final CO₂ and temperature equilibration.

Method 2: The day before embryo culture:
1. Prepare microdrops of supplemented HTF in culture dishes, under oil. Make certain that the oil dishes have had 48 hours to equilibrate prior to making the microdrops and/or use oil from a flask that has been equilibrated with culture medium in the incubator. Using equilibrated oil will ensure that the pH of the medium is reached prior to the onset of embryo culture.

2. Place the culture dishes in the incubator overnight for CO₂ and temperature equilibration.

Preparation of Microdrops for Culture of Embryos from Day 1 to Day 2 or Day 3

Method 1: The day before embryo culture:
1. Label the bottoms of an appropriate number of culture dishes with the patient's name.

2. Using a serological pipette, fill each dish with washed, room-temperature appropriate oil.

3. Place all dishes into the appropriate patient incubator for overnight equilibration.

On the day of embryo culture (Day 1):

1. Use a sterile pipette or tip, to underlay an appropriate number of 20-30 µl microdrops of supplemented HTF into each equilibrated oil dish.

2. Place all dishes into their inverted lids and place in the incubator for 3 to 4 hours.

3. Wash the zygotes and transfer them to the microdrops.

4. Place the dishes in the appropriate patient incubator for culture until Day 2 or Day 3.

5. Evaluate the embryos using your standard laboratory criteria and protocols.

Method 2: The day before embryo culture:

1. Label the bottoms of an appropriate number of culture dishes with the patient's name.

2. Using a serological pipette, fill each dish with washed, room-temperature appropriate oil. Use a sterile pipette or tip, to underlay an appropriate number of 20-30 µl microdrops of supplemented HTF into each equilibrated oil dish.

3. Place all dishes in the incubator overnight equilibration.

On the day of embryo culture (Day 1):

1. Wash the zygotes and transfer them to the microdrops.

2. Place the dishes in the appropriate patient incubator for culture until Day 2 or Day 3.

3. Evaluate the embryos using your standard laboratory criteria and protocols.

SPECIMENS

STERILE A	RX Only	REF	LOT	■	■	■
Sterile Using Asptic Processing Technics	By Prescription Only	Catalogue Number	Batch Code	Consultez le mode d'emploi	Manufacturer	Keep Away From Sunlight
1°C-8°C	EC REP					
Temperature Limitation	Authorized Representative in the European Community	Use By	GST DataMatrix Barcode	Do Not Resterilize	Do Not Use If Package Is Damaged	

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