Instructions for the Use of the 38Special GPS® Dish

(Catalogue Number: SP38-010)

PRECAUTIONS AND WARNINGS

- Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
- Caution: The user should read and understand the Instructions for Use, Precautions and Warnings, and be trained in the correct procedure before using the 38Special GPS® Dish.
- Do not use the product if the product packaging appears damaged or broken.
- For single use only. Do not use after expiry date.
- To avoid problems with contamination, practice aseptic techniques.

GENERAL INFORMATION

Indications for Use

The 38Special GPS® micro-wells dish is specifically designed for oocyte and embryo handling and culture. This 38mm dish replicates the 60 mm dish used in the clinical laboratory and is designed to support the performance of the 38Special GPS® function. The sloped concave micro-well bottoms of the 38Special GPS® dish allow oocytes and embryos to sit at a central location away from the well walls. The dish is a 12 round well pattern in 4 rows. Two large triangle shaped wells are closer to the dishes exterior vertical walls. The triangle wells are by design recommended for washing oocytes and embryos. The concave nature of the wells provides the thinnest well bottom possible, helping to reduce refraction and allow for optimal visualization. The wells may reduce droplet collapsing/mixing, offer better orientation/optics, and reduce set-up/observation time.

Storage and Shelf Life

Store at room temperature. Five (5) years from date of manufacture.

Composition

The 38Special GPS® dish is constructed of polystyrene and is non-pyrogenic. It has passed USP class VI testing and is sterilized by gamma irradiation.

QUALITY CONTROL SPECIFICATIONS

Parameter	Specification
LAL Endotoxin	< 20 EU/device
1-cell Mouse Embryo Assay (% blastocysts at 96 h of culture)	≥ 80%

INSTRUCTIONS FOR USE

The procedures described below have been found to be effective for the handling and culture of human oocytes and embryos in the 38Special GPS® dish. Every laboratory must define and optimize its own procedures.


Off-gassing of the 38Special GPS® dish

38Special GPS® dish packages may be off-gassed with the sleeve partially opened or the dishes may be removed entirely from the packaging. In either case, the opened sleeves, or dishes, must be kept in a sterile environment at ambient temperature and humidity. The dishes should remain covered with the lids until they are loaded with the culture medium and oil.

Marking the 38Special GPS® Dish for Embryo and Patient Identification

Patient identification can be written using a waterproof non-toxic or diamond pen on the outer side of the dish. Alternatively, a label can be placed at the bottom of the dish over some of the drops. (see Figure 1).

Figure 1. Example of the 38Special GPS® dish.

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A. Oil-overlay Technique

1. Carefully pipette 20 µl of culture medium into each of the inner micro-wells, and 60 µl of culture medium into each of the triangle wells. Be careful not to touch the bottoms of the wells with the pipette in order to prevent scratches and/or release of plastic debris into the medium.

2. Pipette 4 ml of oil into the central well covering the micro-wells, using a location away from the wells to gently distribute the oil overlaying the wells previously filled with culture medium.

B. Oil-underlay Technique

1. Pipette 4 ml of oil into the central well covering the micro-wells.

2. Through the oil, carefully pipette 20 µl of culture medium into each of the inner micro-wells and 60 µl of culture medium into each of the triangle wells. Be careful not to touch the bottoms of the wells with the pipette in order to prevent scratches and/or release of plastic debris into the medium.

Pipetting Oocytes and Embryos into the 38Special GPS® Dish

The GPS designed well bottoms have a gently sloping concave bottom which results in the oocytes and embryos migrating towards a central location by gravity. This is generally accomplished within a few minutes as the oocytes and embryos settle after a dish is moved or handled. This feature helps promote rapid visualization of the oocyte or embryo.

Oocytes and embryos can be pipetted in and out of the wells of the 38Special GPS® dish using any type of pipetting device or technique commonly used in the IVF laboratory. Care should be taken not to scratch the well bottoms during pipetting. Generating plastic debris from scratches may adhere to eggs or embryos and scratches may hinder visualization.

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