## English

## PureCepction™

For laboratory procedures only; other uses must be qualified by the end user.		
REF Number	Unit Size	
PureCepction™ 4-determination Kit	ART-2004	3 x 12 mL
PureCepction™ 16-determination Kit	ART-2016	12 x 12 mL

## INTENDED USE

These products are intended to facilitate highly efficient separation of motile sperm from the ejaculate of most semen specimens.

## PRODUCT DESCRIPTION

A 40% (v/v), or 80% (v/v) sterile colloidal suspension of silica particles stabilized with covalently bound hydrophilic silane supplied with HEPES-buffered Human Tubal Fluid (HTF). The components of this kit will allow for the highly efficient separation of motile sperm from the ejaculate of most semen specimens. It is fast, cost-effective, and needs very little equipment or expertise to perform. This product contains 10 mg/L of genta-micin, an aminoglycoside antibiotic.

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## MATERIALS PROVIDED IN THE SPERM SEPARATION MEDIA KIT

- Lower Phase Gradient: 80% PureCepction™ with HEPES-buffered Human Tubal Fluid (HTF-HEPES).
- Upper Phase Gradient: 40% PureCepction™ with HEPES-buffered Human Tubal Fluid (HTF-HEPES).
- Sperm Washing Medium: HEPES-buffered Human Tubal Fluid with Human Serum Albumin, 5 mg/mL.

## PRECAUTIONS AND WARNINGS

PureCepction™ comes packaged with tamper-proof seals and caps. If the seal is broken or the cap loose, do not use product. PureCepction™ has a naturally cloudy appearance. Do not use product if it shows any evidence of particulate matter or contamination. This may be evident by extreme cloudiness or discoloration.

PureCepction™ should remain tightly capped when placed in a CO<sub>2</sub> incubator to avoid pH changes.

**Caution:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found

negative when testing for antibodies to HIV-1/HIV-2, HCV and non-reactive for HBsAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products

prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured by European Pharmacopoeia specifications by established processes.

**Single Use:** To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed. Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

The Sperm Washing Medium (ART-2005-12) component of this kit contains 5 mg/mL human serum albumin.

The products are aseptically processed and supplied sterile.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic gentamicin.

nicin sulphate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

## QUALITY ASSURANCE

One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL. A Certificate of Analysis is available for this product.

## QUALITY CONTROL

All PureCepction™ solutions are tested using one-cell mouse embryo culture and endotoxin assays to ensure quality and safety. However, it is recommended that in-house quality control be performed with each lot.

## USER QUALITY CONTROL

- Examine the PureCepction™ solutions. Do not use any medium that appears discolored. Upper Phase (40%) and Lower Phase (80%) PureCepction™ should have normal opalescent appearance.
- To avoid contamination:
  - Enter all bottles using sterile technique within a suitable sterile environment.
  - Do not use the same sterile pipette or needle to reenter a bottle of medium.
  - When entering a bottle through the stopper via a needle, use a sterile needle. Wipe the stopper with alcohol and air dry. Use a new sterile needle for each gradient.

## MATERIALS REQUIRED BUT NOT INCLUDED

- Sterile disposable polystyrene 15 mL conical centrifuge tube.
- Sterile 5 mL disposable pipettes (graduated) transfer or 3 cc syringes with 1.5"21g needles.
- Centrifuge/level or horizontal. Must be able to rotate for up to 30 minutes at 250 x g to 750 x g.
- 37 °C incubator or water bath.

- Counting chamber.
- Microscope with x10 and x20 objectives.

## DIRECTIONS FOR USE FOR PROCEDURE FOR FRESH SEMEN SAMPLES

- Bring all components of the kit and the semen sample to room temperature before using. This will avoid a cold "shock" to the spermatozoa.
- To prepare the single layered gradient:
  - Transfer 2.0 mL of Lower Phase (PureCepction™ 80%) into the conical tube.
  - Using a pipette or syringe, layer 2.0 mL of Upper Phase (PureCepction™ 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.
  - Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly dispense the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single layered gradient is stable for 1 hour.

**NOTE: Use the PureCepction™ gradient within 1 hour after creating it. Otherwise, the two phases will eventually blend into each other and a sharp interface will not exist. Considerable particulate debris accumulates at the Upper/Lower Phase interface. This is considered a valuable part of the PureCepction™ system.**

- Gently place up to 2.5 mL of liquefied semen onto the Upper Phase using a pipette or syringe.
- Centrifuge for 20 minutes at 350 x g to 400 x g (or up to 750 x g for highly viscous samples). When first centrifugation is completed, you may not be able to see a pellet but it is

- essential to continue the procedure.
- Following centrifugation, remove all layers except the lowest portion (approximately 0.3 mL).
- Add 2 to 3 mL of Sperm Washing Medium and resuspend the pellet.
- Centrifuge for 4 to 8 minutes at 250 x g. The higher number of sperm will require the maximum 8 minute centrifugation time to ensure a complete and thorough sperm wash.
- Remove supernatant and resuspend the pellet in 0.5 mL of Sperm Washing Medium. Sample is now ready for recount and estimate of motility.

## PROCEDURE FOR FROZEN SPERM SAMPLES

Frozen-thawed sperm tend to be less vigorous than fresh sperm. For this reason, use 1.5 mL of Upper and Lower Phases instead of 2.0 mL volumes recommended for fresh sperm.

- Bring all components of the kit and the semen sample to room temperature before using. This will avoid a cold "shock" to the spermatozoa.
- To prepare the single layered gradient:
  - Transfer 1.5 mL of Lower Phase (PureCepction™ 80%) into a conical tube.
  - Using a pipette or syringe, layer 1.5 mL of Upper Phase (PureCepction™ 40%) on top of the Lower Phase. This is done by carefully contacting the surface of the Lower Phase with the tip of the pipette or syringe needle.
  - Allow only the tip of the pipette or bevel of the needle to remain submerged. Slowly dispense the Upper Phase by lifting the pipette or needle up the side of the tube as the level of the Upper Phase rises. An interface separating the two layers will be observed. This single layered gradient is stable for 1 hour.

**NOTE: Use the PureCepction™ gradient within 1 hour after creating it. Otherwise, the two phases will eventually blend into each other and a sharp interface will not exist. Considerable particulate debris accumulates at the 40%/80% interface. This is considered a valuable part of the PureCepction™ system.**

- Gently place the thawed semen sample onto the Upper Phase using a pipette or syringe.
- Centrifuge for 20 minutes at 350 x g to 400 x g. When first centrifugation is completed, you may not be able to see a pellet but it is essential to continue the procedure.
- Following centrifugation, remove all layers except lowest portion (approximately 0.5 mL).
- Add 2 to 3 mL of Sperm Washing Medium and resuspend the pellet.
- Centrifuge for 8 minutes at 250 x g.
- Remove supernatant and resuspend the pellet to 0.5 mL of Sperm Washing Medium.
- Sample is now ready for estimate of count and motility.

**Each laboratory should make its own determination of which medium and protocol to use for each particular procedure.**

Information on specific aspects of IVF and embryo culture is available in our Product Catalog.

## TROUBLESHOOTING

- Occasionally, samples that do not liquefy properly and remain too viscous to pass through the gradient will be encountered. Increasing the centrifugal force up to but no more than 750 x g will aid in separating the sperm in these cases. Pelleting motile sperm from highly viscous semen does not usually present a problem using the PureCepction™ system.

- The most important semen parameter contributing to a pellet of highly motile sperm is sperm progression. The higher the number of progressively motile sperm in the initial sample, the greater the number of sperm you will have in the final pellet. If sperm percent motility or progression is below WHO guidelines for normal semen parameters, you may be able to compensate for shortcomings in the initial semen analysis by using 1.0 mL of both 40% and 80% PureCepction™ instead of 2.0 mL of each. You may also want to leave the last 0.5 to 0.7 mL of the postcentrifugation Lower Phase above the pellet and wash the sperm caught in transit through this portion of the Lower Phase. If this is done, be sure to double the volume of Sperm Washing Medium in order to dilute out this remaining Lower Phase material.

## STORAGE INSTRUCTIONS AND STABILITY

PureCepction™ is assembled in 4- and 16-determination kits. Each determination uses 2.0 mL of Upper Phase (40%) and 2.0 mL of Lower Phase (80%). Sperm Washing Medium may also be supplied.

ART-2004 is assembled with 1 x 12 mL PureCepction™ 40% (v/v), 1 x 12 mL PureCepction™ 80% (v/v), and 1 x 12 mL Sperm Washing Medium.

ART-2016 is assembled with 4 x 12 mL PureCepction™ 40% (v/v), 4 x 12 mL PureCepction™ 80% (v/v), and 4 x 12 mL Sperm Washing Medium.

All PureCepction™ components must be stored at 2 °C to 8 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.

C. Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

## RELATED PRODUCTS


ART-1005/1006 Quinn's™ Sperm Washing Medium  
PureCepction™ gradients are also available in 100 mL volumes (REF # ART-2040, 2080, 2100).  
SAGE™ In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.  
Quinn's™ is registered by CooperSurgical, Inc. PureCepction™ is a trademark of CooperSurgical, Inc.

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