**EC REP** CopperSurgical Distribution B.V.  
Celsiusweg 35, 5928 PR Venlo, The Netherlands  
LifeGlobal Group, LLC, 393 Soundview Rd, Guilford, CT 06437 US  
T: 1800-648-1151 Int'l: +45 4690200  
sales@coopersurgical.com Orders@origo.us

## EN Instructions for the Use of Washing mHTF w/ HEPES & Phenol Red

(Catalogue Numbers: LGFR-100, LGFR-500)

### PRECAUTIONS AND WARNINGS

- Caution:** Federal Law (USA) restricts this device to sale by or on the order of a physician (or properly licensed practitioner).
- Caution:** The user should read and understand the Instructions for Use, Precautions and Warnings, and be trained in the correct procedure before using the Washing mHTF w/ HEPES & Phenol Red when handling gametes and embryos outside a CO<sub>2</sub> incubator.
- Not to be used:** injection.
- Do not resterilize:**
- Do not use the product if:**
  - the product packaging appears damaged or if the seal is broken
  - the expiry date has been exceeded
  - the product becomes discolored, cloudy, or shows evidence of particulate matter
  - the product is not sensitized to the antibiotic
- Washing mHTF w/ HEPES & Phenol Red contains only a low concentration of sodium bicarbonate and should not be gassed, or used for culture under CO<sub>2</sub>.
- To avoid problems with contamination, practice aseptic techniques and discard minimal amounts of excess medium remaining in the bottle.
- Discard unused medium within 7 days of opening.

### GENERAL INFORMATION

#### Indications for Use

For human oocyte retrieval and washing outside of a CO<sub>2</sub> incubator.

#### Storage and Shelf Life

Store at 2-8°C and protected from light. Fourteen (14) weeks from the date of manufacture.

#### Composition

A HEPES-buffered medium is required when gametes and embryos are handled outside of a CO<sub>2</sub> incubator.

Sodium Chloride, Potassium Chloride, Calcium Chloride, Potassium Phosphate, Magnesium Sulfate, Sodium Bicarbonate, Glucose, Lactate Na Salt, Sodium Pyruvate, Phenol Red, HEPES, Gentamicin Sulfate\* (10 µg/ml)

\*from therapeutic-grade source material

### QUALITY CONTROL SPECIFICATIONS

Assay (performed for each batch)	Specification
<b>Physicochemical Tests</b>	
pH	7.2-7.4
Osmolality	280-292 mOsm
<b>Biological Tests</b>	
Endotoxin (LAL)	< 0.5 EU/ml
Sterility Test (bacterial and fungal screen, SAL 10 <sup>-3</sup> )	PASS
<b>Biological Assays</b>	
1-cell Mouse Embryo Assay (% expanded blastocysts at 96 h of culture after 1 h exposure)	> 80%

### INSTRUCTIONS FOR USE

#### Media Required for Oocyte Retrieval

- Needle-Wash Medium: Use Washing mHTF w/ HEPES & Phenol Red to wash the collection needles prior to oocyte retrieval.
- Follicle-Wash Medium: Use Washing mHTF w/ HEPES & Phenol Red to flush the follicles, if appropriate.
- Oocyte-Wash Medium: Use Washing mHTF w/ HEPES & Phenol Red supplemented with protein for the final wash and culture of the oocytes prior to fertilization.
- Oocyte Culture Medium: Use global® for Fertilization supplemented with protein for the final wash and culture of the oocytes prior to fertilization.

The procedures described below have been found to be effective for human oocyte retrieval and washing outside of a CO<sub>2</sub> incubator. Every laboratory must define and optimize its own procedures.

After each time the original bottle is opened recap the bottle tightly and store at 2-8°C, protected from light.

- The day before oocyte retrieval, aliquot appropriate amount of oil into a 14 mL Falcon tube, and place it in a warming block (36-38°C) to equilibrate overnight.
- Twenty-four (24) hours prior to the use of Washing mHTF w/ HEPES & Phenol Red, prepare oocyte-wash medium by supplementing the medium with either Human Serum Albumin (HSA) or LifeGlobal® Protein Supplement to achieve desired % (v/v) of protein supplementation. Store overnight at 2-8°C.

Place a tube of oocyte-wash for each patient in the warming block.

- Place 10-12 100 mL Falcon dishes and 1-60 mL Falcon dishes on the warming surface of the laminar flow hood.
- Just prior to oocyte retrieval, prepare the initial oocyte wash by pouring of oocyte-wash medium into a warmed 60 mm Falcon dishes. Cover the oocyte-wash medium with pre-equilibrated oil and place the dish on the warming surface of the laminar flow hood.

During oocyte retrieval, search for the oocytes under a dissecting microscope. Using a sterile pipette, transfer each oocyte into the initial wash dish containing oocyte-wash medium.

- Dissect each oocyte using a 25 g 1/12" needles to remove the degenerated cumulus cells, debris and blood.

Transfer all oocytes to a dish of the oocyte culture medium for final washing. Pipette all the oocytes a few times to remove traces of oocyte-wash medium.

- At the conclusion of oocyte retrieval, transfer each dissected oocyte into fresh culture media, according to your standard laboratory practice, and place in the incubator until fertilization or ICSI.

### SYMBOLS

STERILE A	RX Only	REF	LOT	■	■	■
Sterile Usando Procedimientos asépticos	By Prescription Only	Catalogue Number	Batch Code	Consult Instructions For Use	Manufacturer	Keep Away From Sunlight

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