HBA™ Tahlili ve PICSİ® Kabi

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1. WHO laboratory manual for the examination and processing of human semen. Sixth edition. 2021
2. Tan J, Taskin O, Albert A, Bedaiwy M. Association between sperm DNA fragmentation and idiopathic recurrent pregnancy loss: a systematic review and meta-analysis. *Reprod Biomed Online* 2019 38(6):951-960
3. Miller D, Pavitt S, Sharma V, Forbes G, Hooper R, Bhattacharya S, Kirkman-Brown J, Coomarasamy A, Lewis S, Cutting R et al. Physiological, hyaluronan-selected intracytoplasmic sperm injection for infertility treatment (HABSelect): a parallel, two-group, randomized trial. *Lancet* 2019; 393:416-422.
4. Jakab A, Sakkas D, Delpiano E, Cayli S, Kovanci E, Ward D, Revelli A, Huszar G. Intracytoplasmic sperm injection: a novel selection method for sperm with normal frequency of chromosomal aneuploidies. *Fertil Steril* 2005; 84(6):1665-73.
5. Huszar G, Jakab A, Sakkas D, Ozenci C, Cayli S, Delpiano E, Ozkavukcu S. Fertility testing and ICSI sperm selection by hyaluronic acid binding: clinical and genetic aspects. *Reprod BioMed Online* 2007; 14(5): 650-663
6. West R, et al. Sperm selection with hyaluronic acid improved live birth outcomes among older couples and was connected to sperm DNA quality, potentially affecting all treatment outcomes. *Hum Reprod* 2022;37: 1106-25
7. Scaruffi P, Bovis F, Casciano I, Maccarini E, De Leo C, Gazzo I, Massarotti C, Sozzi F, Stigliani S, Anserini P. Hyaluronic acid-sperm selection significantly improves the clinical outcome of couples with previous ICSI cycles failure. *Andrology* 2022; 10(4):677-685.
8. Worriolow KC, Eid S, Woodhouse D, Perloe M, Smith S, Witmyer J, Ivani K, Khoury C, Ball G D, Elliot T, Lieberman J. Use of hyaluronan in the selection of sperm for intracytoplasmic sperm injection (ICSI): significant improvement in clinical outcomes—multicenter, double-blinded and randomized controlled trial. *Hum Reprod* 2013;28: 306-14.