**English****Quinn's Advantage™ Thaw Kit**

For laboratory procedures only; other uses must be qualified by the end user.

REF	Product Description	Number	Unit Size
	Quinn's Advantage™ Thaw Kit	ART-8016	3 x 12 mL

**INTENDED USE**

This product is intended for thawing frozen pronuclear and cleavage-stage embryos and blastocysts. This kit is designed to be used in conjunction either with the SAGE™ Quinn's Advantage™ Embryo Freeze Kit (ART-8014) or Quinn's Advantage™ Blastocyst Kit (ART-8015).

**PRODUCT DESCRIPTION**

The components of this kit will allow for the efficient thawing of pronuclear- and cleavage-stage embryos and blastocysts. The components and recommended procedures are the preferred method for improved embryo survivability.

This product contains 10 mg/L of gentami-cin, an aminoglycoside antibiotic.

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**MATERIALS PROVIDED IN THE THAW KIT**

- 1 x 12 mL vial of 0.5 M Sucrose Thawing Medium (REF #8005-12) with 12 mg/mL Human Serum Albumin
- 1 x 12 mL vial of 0.2 M Sucrose Thawing Medium (REF #8007-12) with 12 mg/mL Human Serum Albumin
- 1 x 12 mL vial of Freeze/Thaw Diluent Solution (REF #8013-12) with 12 mg/mL Human Serum Albumin

**PRECAUTIONS AND WARNINGS**

Do not use medium that shows evidence of particulate matter, cloudiness, or is not rose colored.

The 0.5 M Sucrose Thawing Medium (ART-8005-12) component in this kit contains 12 mg/mL human serum albumin.

The 0.2 M Sucrose Thawing Medium (ART-8007-12) component in this kit contains 12 mg/mL human serum albumin.

The Freeze/Thaw Diluent Solution (ART-8013-12) component in this kit contains 12 mg/mL human serum albumin.

To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

**Caution:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when testing for antibodies to HIV-1/ HIV-2, HCV and non-reactive for HBSAg, HCV RNA and HIV-1 RNA. No known test methods can offer assurances that products derived from human blood will not transmit infectious agents. Donors of the source material have been screened for Creutzfeldt-Jacobs disease (CJD). Based on effective donor screening and product manufacturing processes, it carries an extremely remote risk for transmission of viral diseases. A theoretical risk for transmission of

CJD is also considered extremely remote. No cases of transmission of viral diseases or CJD have ever been identified for albumin.

Standard measures to prevent infections resulting from the use of medicinal products prepared from human blood or plasma include selection of donors, screening of individual donations and plasma pools for specific markers of infection and the inclusion of effective manufacturing steps for the inactivation/removal of viruses. Despite this, when medicinal products prepared from human blood or plasma are administered, the possibility of transmitting infective agents cannot be totally excluded. This also applies to unknown or emerging viruses and other pathogens. There are no reports of proven virus transmissions with albumin manufactured to European Pharmacopoeia specifications by established processes.

**Single use:** To avoid problems with contamination, handle using aseptic techniques and discard any excess product that remains in the bottle or vial after procedure is completed.

Reproductive media products are intended for single use only. Re-use of reproductive media may result in using a product past its labeled expiration date or increase the risk of microbial contamination in a subsequent procedure if the practitioner fails to utilize adequate aseptic techniques.

Use of expired or microbial contaminated product may result in suboptimal conditions to promote fertilization and/or embryo quality during in-vitro culture. These conditions may result in the failure of the embryo to develop properly or to implant, potentially leading to a failed assisted reproductive procedure.

**Note:** Embryo is considered a general term. More precisely, SAGE™ considers the period of time initiating when a single diploid cell results from the fusion of male and female

genome resulting in zygote formation with subsequent development from repeated mitotic divisions forming a solid mass or morula (typically day 4-5) and after which a fluid-filled cavity develops resulting in blastocyst formation (typically day 5-6) ending with embryo implantation that begins the end of the first week and is completed by the end of the second week post conception.

**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

This product contains the antibiotic gentamicin sulfate. Appropriate precautions should be taken to ensure that the patient is not sensitized to this antibiotic.

**QUALITY ASSURANCE**  
One-cell MEA tested and passed with 80% or greater blastocyst. USP Endotoxin tested and passed with <1 EU/mL.

A Certificate of Analysis is available for this product.

**DIRECTIONS FOR USE FOR CRYOPRESERVATION OF EMBRYOS**

Controlled hyperstimulation of women undergoing IVF or GIFT produces, on average, 10 to 12 mature oocytes for insemination. It is prudent to replace only a limited number of the resulting embryos, as multiple pregnancies can arise if too many embryos are replaced. Therefore, the majority of patients will have supernumerary embryos. These embryos can be cryopreserved and stored for later use, thus avoiding the necessity of the couple to undertake another stimulated cycle to recover more oocytes for IVF.

The major cause of cell damage during cryopreservation is the formation of intracellular ice during freezing and thawing. By using cryoprotectors, controlling the rates of freezing and thawing, and carefully diluting the cryoprotective from the embryo after thaw-

ing, methods have been developed that allow 80% or more of frozen-thawed embryos to survive and be replaced into the reproductive tract of the woman who produced the oocytes or a genetically unrelated recipient.

**THAWING PROTOCOL**

If the straws have been transferred to liquid nitrogen after being slow-cooled to between -30 °C and -37 °C, they should be thawed rapidly (at least 275 °C/min) so that intracellular ice is swiftly dispersed. This will help prevent cellular damage by the ice crystals. The easiest way to achieve this is to initially hold the straw in air for 30 to 40 seconds and then immerse it in a water bath at 30 °C to 35 °C until the ice has fully melted. This method allows for any liquid nitrogen that may have entered the straw through an imperfectly sealed plug to be blown off before the straw is placed in a water bath. Little loss of straws or their contents occurs with this technique.

Thaw only one cryocontainer at a time. Transfer the liquid contents of the thawed solution to a dry dish and quickly locate the embryos. Pick the embryos up in a minimal amount of solution and transfer them first to 3 mL of 0.5 M Sucrose Thawing Medium (REF #8005-12) at 37 °C for 5 minutes, followed by 3 mL of 0.2 M Sucrose Thawing Medium (REF #8007-12) at 37 °C for 10 minutes, using a new transfer pipette for each procedure to minimize the carry-over of cryoprotectant from one solution to the next. It is recommended that the media be covered with Sterile Oil for Tissue Culture (REF #4008) during use to minimize evaporation of water and a subsequent change in osmolality of the solutions. The embryos are then washed through 7 drops of Freeze/Thaw Diluent Solution (REF #8013-12) at 37 °C. This can be achieved by placing 7 drops, each of 100 µL, under Sterile Oil for Tissue Culture (REF #4008) in a large culture dish. The embryos are placed in each drop and thoroughly washed by pipetting

up and down several times over a period of about 1 minute before being transferred to the next drop. A new transfer pipette should be used after the first drop but the same pipette can be used for subsequent transfers. After the sixth washing drop, the embryos can be transferred to the seventh drop and held for up to 30 minutes at 37 °C before transfer, or placed into culture.

**Each laboratory should make its own determination of which medium to use for each particular procedure.**

Information on specific aspects of IVF embryo culture, and cryopreservation is available in our Product Catalog.

**STORAGE INSTRUCTIONS AND STABILITY**  
Store unopened containers refrigerated at 2 °C to 8 °C. Warm to incubator (37 °C) temperature prior to use. Do not freeze or expose to temperatures greater than 39 °C. The product is stable until the expiration date shown on the label.

- Remove desired volume of product using aseptic procedures.
- Once removed, do not return any volume of product to the original container.
- Do not use if the product becomes discolored, cloudy, turbid, or shows any evidence of microbial contamination.

**RELATED PRODUCTS**

ART-8014 Quinn's Advantage™ Embryo Freeze Kit

ART-8015 Quinn's Advantage™ Blastocyst Kit

SAGE In Vitro Fertilization™ has a full line of products for the Reproductive Medicine Specialist. Please call or write for specific information or to receive a copy of our current catalog. For technical questions, or to reach our Customer Service Department, call the SAGE™ Support Line.

Quinn's Advantage™ is a trademark of CooperSurgical, Inc.

**Call the SAGE™ SUPPORT LINE in the U.S.: (800) 243-2974 International: (203) 601-9818**

**EXPLANATION OF SYMBOLS**

	Catalog Number
	Batch Number
	Use By (year, month, day)
	Do Not Reuse
	Temperature Limitation
	Aseptic Technique Sterilization Membrane Filtered (SAL 10 <sup>-3</sup> )
	<b>ATTENTION:</b> See instructions for use.
	Authorized Representative in the European Community.
	Product conforms to the Medical Device Directive 93/42/EEC
	Manufacturer
<b>RX ONLY</b>	U.S. Federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

**SAGE In Vitro Fertilization**  
*a CooperSurgical Company*

SAGE™ In Vitro Fertilization, Inc.  
a CooperSurgical Company  
95 Corporate Drive  
Trumbull, CT 06611 USA

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Quinn's Advantage™ je ochranná známka společnosti CooperSurgical, Inc.