



# High-quality pipettes for IVF procedures



LISR

## TPC ICSI


- Feel more flexible and bendable
- Longer bevel
- Shorter spike
- Longer and parallel taper



LBC

## TPC Biopsy

Designed to reduce the risk of damaging the cells or embryos during biopsy<sup>1</sup>



LHS

## TPC Holding

Square-shaped, flat and polished surface for maximum grip

Easy to open packaging



# TPC™ Micropipettes

## ICSI Micropipettes

Reference no.	Bend Angle	Tip to Bend Distance (mm)	Bevel Angle	Inner Diameter (µm)	Outer Diameter (µm)	Taper Type
LICR-ID4TA30	30°	0.6	24°	4.5-5.0	5.0-6.0	Parallel Taper
LICR-ID4TL0.6	35°	0.6	24°	4.5-5.0	5.0-6.0	Parallel Taper
LICR-ID4	30°	1	24°	4.5-5.0	5.0-6.0	Parallel Taper
LICR-OD5 (No Spike)	30°	1	20°	4.7-5.3	5.0-6.0	Parallel Taper
LICR-BA21TL0.6 (No Spike)	30°	0.6	21°	5.0-6.0	6.0-7.0	Parallel Taper
LICR-TL0.5	30°	0.5	30°	5.0-6.0	6.0-7.0	Parallel Taper
LICR-TA20	20°	1	30°	5.0-6.0	6.0-7.0	Parallel Taper
LISR	30°	1	30°	5.0-6.0	6.0-7.0	Parallel Taper
LICR-SBS (Short Blunt Spike)	30°	1	30°	5.0-6.0	6.0-7.0	Parallel Taper
LICR-TA35	35°	1	30°	5.0-6.0	6.0-7.0	Parallel Taper

## Holding Micropipettes

Type	Reference no.	Bend Angle	Tip to Bend Distance (mm)	Inner Diameter (µm)	Outer Diameter (µm)
Small	LHC-OD75TA35	35°	1	14-18	75-85
Medium	LHC-OD100TA35	35°	1	14-18	90-110
Medium	LHC-OD100ID20	30°	1	18-23	90-110
Large	LHC-ID15	30°	1	14-18	110-130
Large	LHC-TA35	35°	1	19-27	110-130
Large	LHC-TA20	20°	1	20-27	110-130
Large	LHS	30°	1	22-27	110-125
Large	LHC-ID30	30°	1	28-35	110-130
Large	LHC-ID30TA35	35°	1	28-35	110-130

## Biopsy Micropipettes

Type	Reference no.	Bend Angle	Tip to Bend Distance (mm)	Bevel Angle	Inner Diameter (µm)	Outer Diameter (µm)
Flat and Polished	LBC-OD27BA90	30°	1	-	23-27	27-31
Flat and Polished	LBC-OD30BA90	30°	1	-	25-29	29-33
Flat and Polished	LBC-OD32TA35	35°	1	-	28-32	32-36
Flat and Polished	LBC-OD35BA90	30°	1	-	31-35	35-39
Flat and Polished	LBC-OD35TA35	35°	1	-	31-35	35-39
Flat and Polished	LBC-OD40BA90	30°	1	-	35-39	39-43
Flat and Polished	LBC-OD40TA35	35°	1	-	35-39	39-43
Beveled and Polished	LBSR	30°	1	40°	35-39	29-30

## Hatching and Dissection Micropipettes

Type	Reference no.	Bend Angle	Tip to Bend Distance (mm)	Inner Diameter (µm)	Outer Diameter (µm)
Assisted Hatching	LDS (Flat and Polished)	30°	1	10-12	12-14
Zona Dissection	LCS (Sharp Point)	30°	1	-	-

### Reference

1. Nagy, Z.P. et al., 2019. In vitro fertilization: a textbook of current and emerging methods and devices

© CooperSurgical, Inc. The trademarks used herein are the property of CooperSurgical, Inc. All rights reserved.  
Order No. PIP\_PDS\_0003\_ROW\_US\_V2 · May 15, 2024. All information correct at time of print. Specifications are subject to change without notice or obligation on the part of the manufacturer.