



HANDLING DEVICES

The Stripper[®] and EZ-Range



Innovative handling device options for denudation and handling





Pipettors

Type	Description	Reference no.
EZ-Grip®	EZ-Grip micropipettor	7-72-2800
EZ-Grip	EZ-Grip micropipettor with Comfort Rest	7-72-2802
The Stripper	The Stripper micropipettor	MXL3-STR
The Stripper	The Stripper micropipettor with Comfort Grip	MXL3-STR-CGR
The Stripper CC	Micropipettor for cumulus complex transfer	MXL3-STR-CC
The Stripper PGD	Micropipettor for handling blastomeres	MXL3-STR-PGD
The Stripper PGD	Micropipettor for handling blastomeres with Comfort Grip	MXL3-STR-PGD-CG

Ez-Tip®

	Size of ID (µm)	Multipack (1x vial of 20)	Multipack (5x vial of 20)	Individual Pack (Box of 20)	Primary Use	Pipettor Options/Compatibility
	75	7-72-2075/1	7-72-2075/5	7-72-4075/20	Handling blastomeres for PGT or smaller species oocytes (eg mouse)	Stripper MXL3-STR; MXL3-STR-CGR
	125	7-72-2125/1	7-72-2125/5	7-72-4125/20	Oocyte stripping for ICSI	
	135	7-72-2135/1	7-72-2135/5	7-72-4135/20		
	145	7-72-2145/1	7-72-2145/5	7-72-4145/20		
	155	7-72-2155/1	7-72-2155/5	7-72-4155/20		
	170	7-72-2170/1	7-72-2170/5	7-72-4170/20	Handling embryos	EZ-Grip with Comfort Rest 7-72-2800
	200	7-72-2200/1	7-72-2200/5	7-72-4200/20	Handling blastocysts and expanded blastocysts	EZ-Grip without Comfort Rest 7-72-2802
	290	7-72-2290/1	7-72-2290/5	7-72-4290/20		
	600	7-72-2600/1	7-72-2600/5	7-72-4600/20	Handling larger cell masses and cumulus complexes	

Stripper® Tips

	Size of ID (µm)	Multipack (1x vial of 20)	Individual Pack (Box of 20)	Primary Use	Pipettor Options/Compatibility
	50	MXL3-50	MXL3-IND-50	Handling blastomeres for PGT	Stripper PGD, MXL3-STR-PGD
	75	MXL3-75	MXL3-IND-75	Handling blastomeres for PGT or smaller species oocytes (eg mouse)	
	100	MXL3-100	MXL3-IND-100		
	125	MXL3-125	MXL3-IND-125		
	135	MXL3-135	MXL3-IND-135		
	150	MXL3-150	MXL3-IND-150	Oocyte stripping for ICSI	EZ-Grip with Comfort Rest 7-72-2800
	175	MXL3-175	MXL3-IND-175		
	200	MXL3-200	MXL3-IND-200	Handling blastocysts and expanded blastocysts	EZ-Grip without Comfort Rest 7-72-2802
	275	MXL3-275	MXL3-IND-275		
	600	MXL3-600	MXL3-IND-600	Handling larger cell masses and cumulus complexes	Stripper CC MXL3-STR-CC
1000	MXL3-CC	MXL3-IND-CC			

Accessories

Type	Description	Reference no.
EZ-Grip spares kit	Plunger, spacer and seal replacements	5-72-2000
The Stripper O-ring repair kit	O-rings, spacer	MXL3-ORK
The Stripper Plunger replacement	Plunger replacement kit	MXL3-PLG
The Stripper CC O-ring repair kit	O-rings, spacer	MXL3-ORK-CC
The Stripper PGD O-ring repair kit	O-rings, spacer	MXL3-ORK-PGD
The Stripper PGD Plunger replacement	Plunger replacement kit	MXL3-PLG-PGD
The Comfort Grip	Comfort Grip for The Stripper	MXL3-CGR
EZ-Stand	Stand for pipettors	6-72-2810