

## Safety Data Sheet

Document no: SDS-39

---

**Product:** **MediCult Vitrification Cooling** **Page 1 / 2**

**Date first created: 2010/03/18** **Version 5** **Date revised: 2017/10/12**

### 1. Identification of Substance/Preparation and of the Company

Product name:	MediCult Vitrification Cooling	Manufacturer:	ORIGIO a/s
Cat. No.:	1228		Knardrupvej 2
Product use:	For ultra rapid cooling of oocytes, embryos and blastocysts.		DK-2760 Måløv +45 46 79 02 00

---

### 2. Hazard Identification

The product is classified as non-hazardous according to CLP Regulation (EC) No 1272/2008.

---

### 3. Composition / Information on Ingredients

Component	CAS#	Approx. %
Water		>30
Ethylene Glycol	107-21-1	<20
1,2-Propandiol	57-55-6	<20
Sucrose (Except Equilibration Medium)	50-99-7	<20
Human Serum Albumin		<2
Sodium Chloride	7647-14-5	<1
L-glutamine	56-85-9	<1
HEPES free acid	7365-45-9	<1
Sodium Lactate	113-24-6	<1
Sodium Bicarbonate	144-55-8	<0.1
Potassium Chloride	7447-40-7	<0.1
Calcium Chloride	10043-52-4	<0.1
Potassium Phosphate	7778-77-0	<0.01
Magnesium Phosphate	7487-88-9	<0.01
Sodium Pyruvate	113-24-6	<0.01
Gentamicin sulphate	1405-41-0	<0.01

---

### 4. First-Aid Measures

<b>Inhalation:</b>	Not a likely source of exposure.
<b>Ingestion:</b>	Wash out mouth with water. If swallowed consult a physician.
<b>Skin:</b>	Wash with soap and water after each contact.
<b>Eyes:</b>	Flush with copious amounts of water. Assure adequate flushing by separating the eyelids with fingers. If irritation develops, consult a physician.

---

### 5. Fire Fighting Measures

<b>Fire:</b>	Not considered to be a fire hazard.
<b>Extinguishing media:</b>	Dry chemical, foam, carbon dioxide or water spray.

---

### 6. Accidental Release Measures

Absorb on suitable absorbent, such as paper tissue.

---

### 7. Handling and Storage

Store in original container at 2-8°C, protected from light. Do not freeze prior to use.

---

#### ORIGIO a/s

Knardrupvej 2, 2760 Måløv, Denmark Tel.:+45 4679 0200 Fax:+45 4679 0300  
 Customer Service Tel.: +45 4679 0202 Fax:+45 4679 0302 Web:[www.origio.com](http://www.origio.com) Mail:[origio@origio.com](mailto:origio@origio.com)  
 QI 3.14.01.17, annex 1, v. 05

## Safety Data Sheet

Document no: SDS-39

---

**Product:** **MediCult Vitrification Cooling** **Page 2 / 2**

**Date first created: 2010/03/18** **Version 5** **Date revised: 2017/10/12**

---

### 8. Exposure Controls/Personal Protection

Avoid contact with skin and eyes. Wash thoroughly after handling.

---

### 9. Physiological and Chemical properties

Appearance:	A clear, colourless liquid
Odour:	Odourless
Boiling point:	N/A
Vapour Pressure:	N/A
Specific Gravity:	N/A
Vapour Density:	N/A
Melting Point:	N/A
Solubility H <sub>2</sub> O:	Soluble
Evaporation Rate:	N/A
Flash Point:	None

---

### 10. Stability and Reactivity

Stability:	Stable
Conditions to avoid:	None known
Hazardous Decomposition or by-products:	None known
Hazardous polymerisation:	Not expected to occur

---

### 11. Toxicological Information

No information available.

---

### 12. Ecological Information

No known ecological damage.

---

### 13. Disposal Considerations

Dispose of as medical waste.

---

### 14. Transportation Information

Classified as non-hazardous for transport.

---

### 15. Regulatory Information

The product has been evaluated in accordance with CLP Regulation (EC) No 1272/2008. The product has been classified as non-hazardous.

---

### 16. Other Information

The information contained herein is furnished without warranty of any kind. Users should consider these data only as a supplement to other information gathered by them and must take independent determinations of the suitability and completeness of information from all other sources to assure proper use and disposal of these materials and the safety and health of employees and customers. ORIGIO a/s shall not be held liable for any damage resulting from handling or from contact with its products.

#### ORIGIO a/s

Knardrupvej 2, 2760 Måløv, Denmark Tel.: +45 4679 0200 Fax: +45 4679 0300  
Customer Service Tel.: +45 4679 0202 Fax: +45 4679 0302 Web: [www.origio.com](http://www.origio.com) Mail: [origio@origio.com](mailto:origio@origio.com)  
QI 3.14.01.17, annex 1, v. 05