

The Stripper®

Your trusted brand for
denudation and handling



The Stripper®


The Stripper® Tips

Size of ID (µm)	Multipack (1x vial of 20)	Individual Pack (Box of 20)	Primary Use	Pipettor Options/ Compatibility
50	MXL3-50	MXL3-IND-50	Handling biopsied PGT material	Stripper PGD, MXL3-STR-PGD MXL3-STR-PGD-CG
75	MXL3-75	MXL3-IND-75	Handling smaller species oocytes e.g. mouse	
100	MXL3-100	MXL3-IND-100		
125	MXL3-125	MXL3-IND-125		
135	MXL3-135	MXL3-IND-135	Oocyte denudation for ICSI	
150	MXL3-150	MXL3-IND-150		
175	MXL3-175	MXL3-IND-175	Handling embryos	
200	MXL3-200	MXL3-IND-200	Handling blastocysts and expanded blastocysts	
275	MXL3-275	MXL3-IND-275		
600	MXL3-600	MXL3-IND-600	Handling larger cell masses and cumulus oocyte complexes	
1000	MXL3-CC	MXL3-IND-CC		



The Stripper® Pipettors

Type	Description	Reference no.
The Stripper	The Stripper micropipettor	MXL3-STR
The Stripper	The Stripper micropipettor with Comfort Grip	MXL3-STR-CGR
The Stripper CC	Micropipettor for cumulus oocyte complex transfer	MXL3-STR-CC
The Stripper PGD	Micropipettor for handling blastomeres	MXL3-STR-PGD
The Stripper PGD	Micropipettor for handling blastomeres with Comfort Grip	MXL3-STR-PGD-CG



The Stripper® Accessories and Replacement Kits

Type	Description	Reference no.
The Stripper and The Stripper PGD O-ring replacement kit	O-rings, spacer	MXL3-ORK
The Stripper Plunger wire replacement	Plunger wire replacement kit	MXL3-PLG
The Stripper CC O-ring repair kit	O-rings, spacer	MXL3-ORK-CC
The Stripper PGD Plunger replacement	Plunger wire replacement kit	MXL3-PLG-PGD
The Comfort Grip	Comfort Grip for The Stripper and The Stripper PGD	MXL3-CGR

© CooperSurgical, Inc. The trademarks used herein are the property of CooperSurgical, Inc. All rights reserved • Order No. PIP_PDS_0002_ROW_US_V2 • Version 2 • December 19, 2023. All information correct at time of print. Specifications are subject to change without notice or obligation on the part of the manufacturer.