

Everyone is Unique

Continual cryopreservation monitoring from RI Witness™



Confidence, Efficiency & Trust

CooperSurgical®
Fertility and Genomic Solutions


Cryopreservation with RI Witness

RI Witness provides an unbroken chain of custody throughout every stage of a patient cycle that's in the lab or clinic, including cryopreservation. Whether you use straws or vials, with our e-witnessing system you can separately label and track each and every single embryo or gamete throughout its freeze and thaw processes.

Patients going through IVF can be assured that their embryos are tracked at every step. Linked to their unique RI Witness ID, barcodes on all of a patient's cryopreserved materials are never repeated. Multiple straws or vials from the same patient couple are digitally linked together. Misidentification of donor bank gametes, and embryos or gametes transferred from elsewhere is also safeguarded against.

Barcodes are never duplicated – even between clinics. Each barcode holds a unique identifying code.

Each label's unique barcode is linked to the patient's name, unique ID number, date of birth and the embryo number and date of freezing.



Our Brady® printer range

All our RI Witness-compatible Brady printers generate barcode labels for both straws and vials.



BMP-51

Suitable for smaller cryopreservation volumes, this printer is also ideal if you want to switch between labels for straws and vials easily. No need for separate printers. Supported from v4.0 onwards.



BMP-71

An all-round printer with a wide range of label options, from 100 to 500 labels/roll. This printer is ideal for mid-range cryopreservation volumes, where UDI-barcodes are not required. Supported from v3.0 onwards.



BBP-12

A high-volume printer for clinics freezing high volumes of embryos and gametes, using 1500 labels per roll so you don't need to change rolls constantly. Supported from v4.0 onwards.

Ordering codes for our Brady printer range.

Item	Item numbers	Vol	Brady part number	Unique Barcodes	Label size (mm)	Ribbon composition
------	--------------	-----	-------------------	-----------------	-----------------	--------------------

Printer							
Brady BMP51		6-70-501		BMP51			
Labels	Straws	6-70-520	120	M-86-461	Yes	15,24x66,68	Integrated in cartridge
	Vials	6-70-521	120	M-75-46	Yes	25,40x66,68	Integrated in cartridge

Printer							
Brady BMP71		6-70-701		BMP71			
Labels	Straws	6-70-720	250	M71-19-427	Yes	25,40x25,40	6-70-710
		Vials	6-70-721	100	M71-86-461	No	15,24x66,68
	6-70-722		500	M71-76-461	No	15,24x41,28	6-70-711
	6-70-724		250	M71-19-492	Yes	25,40x25,40	6-70-712
	6-70-727		100	M71-115-461	Yes	25,40x66,68	6-70-711
	6-70-728	250	M71-106-461	Yes	25,40x44,45	6-70-711	
Ribbon		6-70-710		M71-R4300		50,8x46	
		6-70-711		M71-R6200		50,8x46	
		6-70-712		M71-R6400		50,8x46	

Printer							
Brady BB12		6-70-901		BBP12			
Labels	Straws	6-70-920	1500	THT-153-461-1.5-SC	Yes	15,24x66,68	6-70-910
	Vials	6-70-921	1500	THT-131-461-1.5-SC	Yes	25,40x66,68	6-70-910
Ribbon		6-70-910		R7950 110mmx70m /O		110x70	

Copyright CooperSurgical, Inc. • Order No. WIT_FLY_010_V2_ROW - October 7, 2020. All information correct at time of print.