

What do people say about RI Witness?



"Since I started using RI Witness in my laboratory the embryologists sleep better and so do the patients."

Jean-Claude Jacquet, Oriade PMA Grenoble, France

"RI Witness has allowed our lab staff to perform procedures with total confidence and without continual interruptions, giving all team members peace of mind."

Carli W Chapman, The Reproductive Medicine Institute
Oak Brook, USA

Why RI Witness?

We know that embryologists are some of the most dedicated, exacting professionals you can meet. They pride themselves in the care and attention they give each and every patient at their clinic.

Since the first use of RI Witness by a clinic in 2007, we have continued to work closely with leaders in the field to constantly improve this support and safeguarding system which helps and guides embryologists throughout the IVF cycle.

It ensures that every step of the fertility cycle, within the lab, is managed properly and securely, maintaining the safety of your samples.

RI Witness is increasingly becoming the gold standard for security in clinics on every continent around the world.



For more information including videos, customer testimonials and clinic locations please visit:

fertility.coopersurgical.com

CooperSurgical Inc.
fertility.coopersurgical.com

CooperSurgical[®]
Fertility and Genomic Solutions

©CooperSurgical • Order No. EQU_FLY_008_V2_ROW • November 7, 2019 • All information correct at time of print. Specifications are subject to change without notice or obligation on the part of the manufacturer.

RI Witness[™]

Safeguarding your eggs, sperm & embryos throughout the IVF cycle



CooperSurgical[®]
Fertility and Genomic Solutions



We understand you want to know everything is being taken care of as you embark on a life changing journey.

What is RI Witness?

RI Witness has been created to help minimize errors that may occur during ART treatment cycles. This electronic witnessing system supports and guides embryologists in their daily work by monitoring sample movement from one dish or tube to another within the laboratory twenty-four-seven.

This system is automatic, so it does not require additional steps by embryologists to identify samples. This means the embryologist can give their full attention to caring for your samples.


What does this mean to you?

Choosing a clinic with RI Witness installed allows you to rest assured that your samples are being managed and tracked in the securest possible way. In fact, the clinic has gone beyond national legislation requirements by investing in a system providing the highest possible security for the safeguarding of your samples.

A clinic with RI Witness places you at the forefront of its priorities and proves that taking extra measures and every precaution is the centre of their philosophy.




How does it work?




1
Each patient is assigned an ID card with their unique personal code

ID tags containing this unique personal code are fixed to all plasticware holding patient samples




4
RI Witness confirms the identity of samples and allows the embryologists to proceed to insemination



3
Every work area in the lab detects wireless signals from these ID tags.

They are identified, tracked and recorded at each step of the ART process



5
In the unlikely event that samples are mismatched, the system alerts the embryologist, stops proceedings and helps prevent any potential errors

6
Prior to embryo transfer the patient's ID card is automatically checked against the tagged embryo dish



RI Witness helps safeguard your samples during ART processes - from start to finish